


Product Name :

Astra Scientific Three Gang Moulds 40x40x160 mm

Product Code :

ASTCEMENT-0031

**Description :**

Astra Scientific Three Gang Moulds 40x40x160 mm

Technical Specification :**General Descriptions:**

Used to produce the mortar specimens for compression and flexural tests. Made from special alloy steel, they are proposed in two version:

- Standard having a minimum surface hardness of HV200 and
- Heavy duty having a minimum surface hardness of HV400 which is recommended by the EN standards. Furthermore, this model is surface protected by a rust inhibitor treatment and is supplied complete with traceable certificate of conformity and identified by a serial number.

Both models can be fit with the type B measuring pegs for determining the coefficient of thermal expansion conforming to EN 1770. See

accessories.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265