


**Product Name :**  
Astra Scientific 50 mm Three Gang Cube Mould

**Product Code :**  
ASTCEMENT-0032



**Description :**

Astra Scientific 50 mm Three Gang Cube Mould

**Technical Specification :**

**General Descriptions:**

The 50 mm cube specimens are produced for compression test. The dimension is compatible with ASTM C109. Made from steel, precisely machined.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

