


Product Name :
Astra Scientific 50 mm Three Gang Cube Mould

Product Code :
ASTCEMENT-0032



Description :

Astra Scientific 50 mm Three Gang Cube Mould

Technical Specification :

General Descriptions:

The 50 mm cube specimens are produced for compression test. The dimension is compatible with ASTM C109. Made from steel, precisely machined.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

