

Product Name :
Cement Curing Cabinet

Product Code :
ASTCEMENT-0034



Description :

Cement Curing Cabinet

Technical Specification :

General Descriptions:

Ideal for curing cement specimens in commercial and site laboratories, very practical and easy to use. The frame is made of strong polypropylene structure, which is chemical resistant and particularly suitable for cement applications. The front doors are fitted with transparent glass. The humidity is maintained from 95% to saturation by water nebulizer. The temperature is maintained to $20\pm 1^{\circ}\text{C}$ by a immersion heater and separated refrigerator unit . The four stainless steel racks of the internal frame can support the moulds with specimens and a large number of cement prisms. It can also be used to concrete cubes and other mortar specimens. Using a simple plastic pan is also possible to cure specimen in water.

The unit has to be connected to a suitable air compressor.

- Power: 1700 W approx.
- Internal dimensions (lxdxh): 1115x435x1500 mm
- External dimensions: 1160x550x1900 mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265