


Product Name :
Astra Scientific Devel Attrition Tester

Product Code :
ASTAGGRE-0007



Description :

Astra Scientific Devel Attrition Tester

Technical Specification :

Specification: It consists of two hollow cylinders closed at one end and provided with Fitting covers at the either end. These cylinders are mounted on a shaft at angle of 300 with the axis of rotation of the shaft. The shaft rotates at 30-33 RPM. Through a reduction gear operated by a motor and is provided with a revolution counter. Complete with Abrasive Charge consisting of 12 Nos. Hardened steel Balls of 48mm dia. Suitable for operation on 440 Volts, Three Phase, 50 Cycles, A.C. supply.

Note : Option of Digital Preset Counter can be provided at an Extra Cost.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com