

**Product Name :**  
Astra Scientific Wide Wheel Abrasion Testing Machine

**Product Code :**  
ASTAGGRE-0008



**Description :**

Astra Scientific Wide Wheel Abrasion Testing Machine

**Technical Specification :**

Specification: Abrasion Testing Machine is designed for determining the resistance to abrasion/wear of natural stones and concrete products used for paving. The abrasion wheel is 70 mm thick and rotates at a speed of 75 rpm. The machine is equipped with a digital counter which stops the machine at the end of a preset number of revolutions.

Narrow Wheel Abrasion Testing Machine is designed for determining the resistance to abrasion of unglazed tiles, grouts used for tiles and clay pavers. The abrasion wheel is 10 mm thick. Abrasive corundum sand, calibration marble (boloneisse) and mould (for EN 12808-2) should be ordered separately. Abrasion Testing Machines can be convert the other.

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



**Astra Scientific**