

Product Name :
Astra Scientific Efflux Index Apparatus

Product Code :
ASTAGGRE-0032



Description :

Astra Scientific Efflux Index Apparatus

Technical Specification :

Specification: Used to measure the efflux index of fine aggregates (shape and angularity). The efflux index of an aggregate is the required time in seconds of a known volume of aggregates to flow from a known opening. The unit is basically formed by aluminium body, 90 mm diameter by 125 mm height, aluminium feed hopper 100 mm dia. by 170 mm height, control shutter, polycarbonate funnel having 85 mm height, 60° conical part, which end has a diameter of 12 mm, base support, valve, decanter.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265