

**Product Name :**  
Astra Scientific Marshall Stability Test Apparatus

**Product Code :**  
ASTBATE-0006



**Description :**

Astra Scientific Marshall Stability Test Apparatus

**Technical Specification :**

Specification: Marshall Apparatus Consists Of:

- Robust construction
- 50KN capacity
- Rate of travel 50.8mm/min
- Safety cut-off switch

It consists of body housing, a geared screw jack and motor drive mechanism suitable for operation on 230 V, 50Hz, Single Phase, A.C. supply.

MARSHALL LOAD FRAME 1 No.

Breaking Head Stability Mould, with a dial gauge ( having 25mm travel and 0.01mm least count), for flow measurement 1 No.

Compaction Mould, Steel cylindrical 3 Nos.

Base Plate 3 Nos.

Extension Collar 3 Nos.

Compaction Pedestal, Manual Operation, Comprising a steel plate capped on a wooden post. A Mould Clamp is fitted to the top of the plate . 1 No.

Compaction Hammer, satisfies BS 598. The hammer has a 4535 g sliding weight with a free fall of 457 mm 2 Nos.

Load Transfer Bar 1 No.

Paper Discs, Non absorbent, 11 cm diameter pack of 100 No. 1 Pkt.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265