

Product Name :
Astra Scientific Centrifuge Extractor Motorised

Product Code :
ASTBATE-0010



Description :

Astra Scientific Centrifuge Extractor Motorised

Technical Specification :

Specification: Centrifuge Extractor, Electrical Operation, Capacity 1500g, with a Dimmer stat for speed control from 2,400 to 3,600 rpm. Suitable for operation on 230 V, 50 Hz, Single Phase, A.C. supply.

Optional Extra: Filter Paper Discs, Set of 25



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265