

Product Name :
Astra Scientific Mastic Asphalt Tester

Product Code :
ASTBATE-0030



Description :

Astra Scientific Mastic Asphalt Tester

Technical Specification :

Specification: For determining the hardness number of Mastic Asphalt for flooring. It consist of an internally insulated cabinet to the base of which is fitted a water bath having two taps. The bath is heated by an immersion heater and the temperature is controlled with the thermostat at $35^{\circ} \pm 0.5^{\circ} \text{C}$. A 6.5mm dia pin is loaded on the specimen by means of a lever to give a 31.7 Kg weight. A dial gauge 0.01 x 25 mm is provided to record the penetration of pin into the specimen.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265