

Product Name :
Astra Scientific Saybolt Two Tube Digital Viscometer

Product Code :
ASTBATE-0045



Description :

Astra Scientific Saybolt Two Tube Digital Viscometer

Technical Specification :

Specification: The Saybolt Viscometer is used to determine empirical measurement of Saybolt Viscosity of petroleum products at specified temperatures.

The viscometer can be used for temperatures between 21 to 99°C (70 to 210°F) the viscometer includes wateroil bath, stirrer, cooling coil, electric heater with digital thermo regulator, furol orifice, universal orifice, thermometer support and 2 x 60 ml glass Saybolt Viscosity Flask.

Viscosity thermometer set consists of 6 thermometers with the temperature ranges; 19 to 27°C, 34 to 42°C, 49 to 57°C, 57 to 65°C, 79 to 87°C (250 mm length) and 95 to 103°C where each thermometer with 0.1°C subdivisions.

Filter funnel, withdrawal tube and thermometer set should

be ordered separately.

Supplied complete with:

Universal orifice

Furol orifice

Thermometer support

Heat transfer oil, 5 lt

Key Saybolt Viscosity Flask, Glass, 60 ml, 2 pcs.

Saybolt Two-tube Digital Viscometer 220-240 v 50-60 Hz

Filter Funnel with Wire Mesh and Clip

Withdrawal Tube

Saybolt Viscosity Thermometer set, 6 Pcs.

Saybolt Viscosity Flask, Glass, 60 ml

Dimensions : 450 x 300 x 550 mm

Weight (approx.) : 10 Kg

Power : 750 W



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265