

Product Name :
Astra Scientific Say Bolt Viscometer

Product Code :
ASTBATE-0046



Description :

Astra Scientific Say Bolt Viscometer

Technical Specification :

Specification: Say bolt Viscometer, Electrically Heated, ASTM D88, D244, AASHTO T72 for the empirical measurement of Say bolt Viscosity of petroleum products at specified temperatures between 70OF and 210OF. This is also used for determining the Saybolt Furol Viscosity of bituminous materials at temperatures of 250, 275, 300, 350, 400 and 450OF. It comprises one each of Cylindrical Oil Cup, Universal Tip, Furol Tip, Bath Fitted with immersion Heater mounted on a stand. Dimmer stat for temperature control, Stirrer with shield. Complete with insulated handle and thermometer support receiving flask, withdrawal tube, filter funnel, thermometer support for cup and circular spirit level. Suitable for operation on 230 V 50 Hz, Single Phase, A.C.

Brass Oil Tube

Receiving Flask, 60ml

Universal and Furol Tip (Set of 2)

Filter Funnel

Thermometer Support

Heating Coil

ASTM Thermometer Type 18 F

Range : 66°F to 80°F

ASTM Thermometer Type 19 F

Range: 94°F to 80°F

ASTM Thermometer Type 19 F

Range: 120°F to 148°

ASTM Thermometer Type 20F

Range: 134°F to 148°F

ASTM Thermometer Type 21 F

Range : 174°F to 188°F

ASTM Thermometer Type 22 F

Range : 204°F to 218°F

ASTM Thermometer Type 77F

Range: 245°F to 265°F

ASTM Thermometer Type 77F

Range: 295°F to 315°F

ASTM Thermometer Type 78 F

Range: 345°F to 365°F

ZI 4019-P ASTM Thermometer Type 80F

Range : 395°F to 415°F

ZI 4019-Q ASTM Thermometer Type 81F

Range: 445°F to 465°F



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265