

Product Name :
Astra Scientific Bar Cutting Machine

Product Code :
ASTSTEEL-0010



Description :

Astra Scientific Bar Cutting Machine

Technical Specification :

TMT Cutting Machines, also known as Bar Cutting Machines are lightweight, easy in operation, yet very sturdy in structure. The dense building of this rebar equipment makes it ideal for a standard construction site.

Durable fast and easy to use, this machine gives a one-shot cutting operation with perfect cutting accuracy thus minimizing the use of labor and long drawn cutting process for single cuts. The extra characteristics of the Rebar Cutting Machine help save a great heap of money and time for the customer in vast quantity, thus increasing their output capability.

This workhorse knows no stop as it can be effectively put to use for a constant round-the-clock production.

Technical Specifications:

R.P.M

Machine Motor

H.P.	Kg/mm245			Kg/mm265			Kg/mm280			Dimensions LxWxH			Weight approx.
	1Ø	2Ø	3Ø	1Ø	2Ø	3Ø	1Ø	2Ø	3Ø				
90	3	28	16	8	25	12	6	20	10	10	835 x 560 x 825 MM	450 KG	
90	5	36	20	10	32	16	8	25	12	6	1050 x 600 x 920 MM	480 KG	
90	5	42	25	12	36	20	10	32	16	8	1200 x 750 x 910 MM	680 KG	
47	7.5	52	32	20	42	25	16	36	20	12	1300 x 710 x 940 MM	770 KG	
40	7.5	55	36	28	48	32	28	44	28	20	1325 x 720 x 940 MM	800 KG	



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265