

Product Name :
Astra Scientific Bar Bending Machine

Product Code :
ASTSTEEL-0011



Description :

Astra Scientific Bar Bending Machine

Technical Specification :

Descriptions:

This machine effectively increases the production capacity of your steel yard, minimizing the use of manual labor. TMT bending machine, also known as, Bar Bending Machine are boisterously designed for life long operation making them uncomplicated to use and extremely user friendly. Our Rebar Bending Machines are sternly designed to reduce the role of manpower as well as maintain low cost of production for our clients. Rebar benders are manufactured in large batches, under tight quality controlled environment and are thoroughly tested before shipping them across the globe.

The automatic angle selection permits precise bend at a preset angle making it one step bending process for various forms of bends and stirrups. Bar Bending Machine is easy to use, the machine can be operated by a layman with minimal experience. Rebar Bending machine is highly appreciated by our clients for its high strength, outstanding performance and for being cost effective. Our Research and Development

department is highly modified with the employment of the latest technology. Therefore, we incline to make changes in the machine to prepare it for easier usage keeping cost under the bindings.

Technical Specifications:

R.P.M Machine	Motor H.P.	Kg/mm245			Kg/mm265			Kg/mm280			Dimensions LxWxH	Weight approx.
1Ø	2Ø	3Ø	1Ø	2Ø	3Ø	1Ø	2Ø	3Ø				
9	5	36	20	10	32	16	8	25	12	6	860 x 910 x 880 MM	500 KG
8	5	42	25	12	36	20	10	32	16	8	960 x 970 x 900 MM	700 KG
6	7.5	52	32	20	42	25	16	36	20	12	960 x 970 x 900 MM	780 KG
5.6	7.5	55	36	28	48	32	28	44	28	20	1100 x 1037 x 892 MM	820 KG



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265