

**Product Name :**  
Astra Scientific Stirrup Bending Machine

**Product Code :**  
ASTSTEEL-0014



**Description :**

Astra Scientific Stirrup Bending Machine

**Technical Specification :**

**Descriptions:**

The stirrup bending machine SLB-14 from is a multipurpose bender. Manufactured with high quality raw materials, these run for a long span of time saving time of repair and cost of maintenance.

The automatic stirrup bending machines from are lubricated properly, which ensures lower rate of friction. This is another reason why the automatic rebar stirrup bending machine from ensures high productivity for a longer time period.

The expert engineers from do not compromise on safety and hence the TMT stirrup bending machines from this leading construction equipment manufacturer are compatible with the aggressive working environment.

We have gained expertise in providing our clients with Stirrup Bending Machine that is available in different models and sizes at most competitive prices. All our machines are manufactured using superior raw

material and latest technology to assure that final output meets with international quality norms & standards. We also check these machines on various parameters to ensure flawlessness.

Technical Specifications:

Length 850 mm

Width 685 mm

Height 885 mm

Weight (approx.) 400 KG

Machine 20 RPM

Max. Bending DIA 12 mm - FE 500 TMT

Power Rating 1.1 kw

Control 415 V, 3 DIA, A.C

Minimum Stirrup Size 100 x 100 mm



## Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265