

Product Name :
Astra Scientific Polishing and Lapping Machine

Product Code :
ASTROCK-0007



Description :

Astra Scientific Polishing and Lapping Machine

Technical Specification :

Specification:

This unit is provided with a 1 HP Single Phase, A.C. motor. This bench-mounted single-spindle lapping machine is ideally suited for the final polishing of mounted rock or concrete specimens. Two sample holders of sample size ordered for are provided with each machine as standard supply. Sample holders of other size can also be provided on request. This is a motor-driven unit with 450/500 rpm. A swing-in tap, for continuous water supply during operation, is also provided.

This unit is operated on Single Phase. 230 V A.C. motor.

- A Sample Holder for Ex size
- B Sample Holder for Ax size
- C Sample Holder for 38mm size
- D Sample Holder for Bx size

E Sample Holder for Nx size

- Compact table model
- 20 cm dia top and Adapter to hold polishing cloth or paste
- Sample Holders to accommodate upto NX size Cores
- Continuous water feed arrangement during operation



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265