


Product Name :
Astra Scientific Digital Point Load Test Apparatus

Product Code :
ASTROCK-0010



Description :

Astra Scientific Digital Point Load Test Apparatus

Technical Specification :

Specification: Digital Point Load Test Apparatus consists of a 60 kN capacity load frame with a hydraulic loading ram driven by a hand pump. The frame is adjustable for testing of samples up to 102mm diameter. A ruler assembled on the frame allows the direct measurement of the distance between the conical platens before and after the test. The compression load is measured by a pressure transducer connected to an advanced digital display unit assuring the best accuracy and resistance to the failure shocks.

Safety Features:

- Load range : 0 - 60 kn
- Digital display : 2 x 16 characters
- Resolution : 32.000 div.
- Accuracy : $\pm 1\%$
- Load pacer included

- Load measurement in both kn and mpa
- Serial port for pc connection



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265