

Product Name :
Astra Scientific Boring Machine

Product Code :
CNCWORK0009



Description :

Astra Scientific Boring Machine

Technical Specification :

The boring machine is mainly used for reboring engine cylinders of motor cycles, automobiles and middle & small-tractors.

--- Reliable performance ,widely use, processing accuracy high productivity.

--- Easy and Flexible operation .

--- Good rigidity ,amount of cutting

Boring diameter 65 – 140 mm 80 – 165 mm

Max. boring depth 300 mm 400 mm

Spindle speed 250, 380 r/min 200, 300 r/min

Spindle feed 0.11 mm/r 0.12 mm/r

Spindle quick reset 636 mm/sec 780 mm/sec

Overall dimensions 320x330x850 mm 360x345x1100 mm

Machine weight 62 kg 85 kg



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265