

Product Name :
Astra Scientific Slotting Machine

Product Code :
CNCWORK0014



Description :

Astra Scientific Slotting Machine

Technical Specification :

This model can produce slotting strokes from 10 to 150mm. The longitudinal and cross movement for ranges about 200 mm and 110 mm respectively. The Ram adjustment is 150mm and the column to Ram distance is 210mm similarly table to Ram distance is 310 mm. This equipment is made of rigid construction and it creates no vibration while machining. It is made of fine grained and seasoned casting that makes the instrument work for long term application.

- Rigid Construction hence smooth and Vibration free operation
- Fine grained castings and duly Seasoned.
- All gears are Machine cut and Properly fitted.

ADJUSTABLE Stroke (mm) 10 to 150 10 to 250

Longitudinal Movement (mm) 200 230

Cross Movement (mm) 110 230

Speed adjustment (Speed) 3 3

Ram Adjustment (mm) 150 250

Motor (HP) 1 1.5

Column to RAM Dist.(mm) 210 330

Table to RAM Dist. (mm) 310 410

Table size (mm) 280x335 280x335

Machine Dimension 1050x550x1175 1150x600x1375

Machine weight (kg) 375 450



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265