

Product Name :

Astra Scientific Block Making Machine Hand Operated

Product Code :

CNCWORK0021

**Description :**

Astra Scientific Block Making Machine Hand Operated

Technical Specification :

Manual Hollow Block Machine is known for their optimum performance and easy operability. These are manufactured using the highest grade of material and are in conformation with the international quality standards. Our manual hollow block machines can also be customized as per our clients specifications and are also available in standard specification.

Features:

High performance

Efficient

Ease of use

Bottom Vibrator 1.5 HP Single / Three Phase Motor

Top Vibrator

Vibration 6000 hz Amplitude Adjustable

Moldings Area 375 X 600 mm

Weight 450 Kg

Pavers Block 4 Block Stroke 800 -1000 block Shift

Hollow Block 100 x 200 x 400 4 Block /Stroke 500-800 Block/Shift

150 x 200 x 400 2 Block /Stroke 250-400 Block/Shift

200 x 200 x 400 2 Block /Stroke 250-400 Block/Shift



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265