

Product Name :

Astra Scientific Fuel Injection and Engine Management Test Bench

Product Code :

ThermalLab0008

**Description :**

Astra Scientific Fuel Injection and Engine Management Test Bench

Technical Specification :

The EFI Test Bench is designed to study the fundamentals of Engine Management System. Different Sensors like speed sensor, MAP, Throttle Position Sensor, Temperature sensor, CAM Sensor etc are used to generate the signal which is sent to the ECU. The EFI system having Electronic Control Unit which actuate the Injector, Ignition system, Duration of Injection, Start of Injection, Ignition timing, Dwell angle etc can be varied. All these signals can be monitored on the oscilloscope. K-Line communication is used to interface between computer and the ECU. The Application software is provided along with the system to monitor all the Input-Output parameters on the computer. Different LUT like DOI, SOI, SOS, Dwell time, Correction factor for Dwell Time and DOI can be optimized by the user under different speed conditions.

Product Specification

Automation Grade Automatic

Frequency (Hz) 50

Speed (rpm) 50-400 rpm

Material Mild Steel

Power (kilowatt) 11 kw

Output Voltage 220V/380V



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265