

Product Name :
Astra Scientific Corollis Component of Acceleration Apparatus

Product Code :
DynamicsLab0004



Description :

Astra Scientific Corollis Component of Acceleration Apparatus

Technical Specification :

The experimental set-up is designed to study Coriolli's Component of Acceleration of a slider crank mechanism. Here the mechanical slider system is replaced by a continuous stream of water flowing through a steady rotating pair of tubes. These tubes can be rotated at various speeds by using a swinging field motor which also acts as dynamometer. A Perspex window on top of the tank gives clear view of the process and prevents splash of the water over the side of the tank. The dynamometer continuously measures torque applied to the rotating tubes.

FEATURES:

- Coriolis components of Acceleration can be determined at various speeds of rotation at water flow rates.
- The equipment's is self contained, water re-circulating, provided with its own speed control unit and water-circulating pump.

SPECIFICATION:

- Water Tank: Made of Stainless Steel
- Rotating Arms: 9mm/6mm orifice diameter, Length 300 mm.
- Rota meter: 250 to 2500 LPH
- Electric Motor: Swinging field, Variable speed.
- Pump: FHP for water circulation.
- The whole set-up is ingeniously designed and schematically arranged on a powder-coated rigid structure
- Required for operation

Water supply and drain

230V, 50Hz, 1 phase

230V, 60Hz, 1 phase; 120V, 60Hz, 1 phase



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265