


**Product Name :**  
Astra Scientific Rack and Pinion Gear Trainer

**Product Code :**  
DynamicsLab0010



**Description :**

Astra Scientific Rack and Pinion Gear Trainer

**Technical Specification :**

The cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections. Many parts have been chromium-plated and galvanized for a longer life.

The model will be mounted on to a tubular frame; Suitable color imported painting will be carried out along with Miracle coating for the model will be done for extra gloss and shine. The painting will be carried out in such a way that different colors will be use for different components such as identification of sectioned area etc. according to the colors code for easy identification of different systems and mechanisms. All the hardware's and gears will be suitably electroplated.

**FEATURES:**

- Operated: Manually (Through a Crank Handle)
- Equipment Layout: On Base
- Fine finish

- Simple operations
- Functional accuracy

SPECIFICATION:

- The properly constructed model is hand Operated,
- Sturdiness
- Corrosion-resistance
- Enhanced durability



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265