

**Product Name :**  
Astra Scientific Cut Section Model of Open Hermetic Compressor

**Product Code :**  
RefrigerationLab0016



**Description :**

Astra Scientific Cut Section Model of Open Hermetic Compressor

**Technical Specification :**

The Cut Section Model Of Open Hermetic Compressor model will be made out of original parts. The model will be suitably sectioned to demonstrate the internal construction details showing the minute information, and working of the same, the model will be suitably painted and mounted on a metal or Mild Steel Sheet base & supplied with key card & very interesting literature regarding working of the same.

**FEATURES:**

- This model will be made out of original parts.
- The model will be suitably painted and mounted on a metal or wooden base

**SPECIFICATION:**

- Kit of a piston compressor
- Compressor consisting of piston and cylinder, housing, driving gear, cylinder cover with valves



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265