

**Product Name :**  
Astra Scientific Freezing Point Apparatus

**Product Code :**  
FuelTesting0030



**Description :**

Astra Scientific Freezing Point Apparatus

**Technical Specification :**

Standard: ASTM D1177

Determines the freezing point of aqueous engine coolant solutions by cooling a sample with continuous agitation until a plateau is observed in a time-temperature curve. Conforms to ASTM D1177 specifications.

**Product Description** The electric freezing point apparatus is a manual freezing point apparatus.

Double tube Ø 30 x 240 mm fitted with cap

- Dewar jar Ø 75 x 280 mm
- Geared motor for stirring at 80rpm with wire stirrer
- Digital thermometer
- Mounted on a plate structure painted with anti-acid epoxidic products

**ACCESSORIES:**

ASTM-114C Thermometer ASTM 114C IP 14C

102-376 Type A Nitrogen collar made in glass

SPARE PARTS:

102-371 Double tube 30 x 240mm

102-372 Wire stirrer

102-373 Dewar jar

102-374 Cork series

102-375 Type B Modified moisture proof collar made in glass



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265