

Product Name :
Astra Scientific Deep Seeker

Product Code :
Civil0001



Description :

Astra Scientific Deep Seeker

Technical Specification :

Deep Search Metal & Mine Detector

1. Light weight, Rugged and professional design
2. Reliable high-speed pulse-induction technology
3. Independent controls for Power, Sensitivity, Threshold and Audio
4. Adjustable parameters allow use as a professional mine or metal detector
5. Fully static detection target metal close to large metals
6. Capable of detecting target metal close to large metals
7. Detects all metals and alloys to a depth of over 1.5 meters with standard coil
8. Continuously variable pitch to indicate size and proximity of metal
9. Bar-graph display for visually in pointing metal

- 10. Concentric multi purpose 3-coil probe mounted on extendable telescopic pole
- 11. Large size coil and stick probes available as accessories for special requirements
- 12. Coils and probes are submersible in water
- 13. Ni/MH, SMF or Dry battery power source with 90-270V internal battery charger



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265