

Product Name :
Astra Scientific PCB Curing Oven Machine

Product Code :
PCB0001



Description :

Astra Scientific PCB Curing Oven Machine

Technical Specification :

PCB Curing Oven Machine For Electrical Lab Training.

PCB Curing Oven Machine
Tabletop unit for curing of liquid protoresist or drying of PCB's.
Timer controlled system allows fast and efficient PCB curing.
PCB size 10" x 12"
Timers : 0-30 minutes.
Temperature : 0-90 degree C



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265

