

Product Name : Insize 150mm Digital Vernier Caliper	Product Code : AST-MES0001
---	--------------------------------------



Description :
Insize 150mm Digital Vernier Caliper
Technical Specification :
Insize is a well-established brand that has launched a digital Vernier caliper with a size of 150 mm. It comes with a manufacturing warranty period of 1 year. It is a basic model that supports data output. With DIN862 specification and 460 g weight, it has emerged as an exclusive product. It supports a power supply of SR 44 or LR 44. The product

dimensions are 11x5x2 which is a standard size. The lower jaw length is 40 mm and it comes with an LCD display for a clear view. It supports a measurement range of 0 to 150 mm. It comes packed in a plastic box. It is made with hardened stainless-steel material. It comes with a zero, power on, and off button. In terms of durability, it is packed with an extended lifespan. It is a great product that is built with quality-assured materials to ensure superior performance. It is a versatile product with stable micro-process slides.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265