

Product Name :
Adapter Flash Chromatography Reservoirs

Product Code :
GLASSWARESUK-AS07LAB004



Description :

Adapter Flash Chromatography Reservoirs

Technical Specification :

Made of borosilicate glass.

Capacity of 250ml.

Chromatography reservoir.

With 24/40 joint.

Made by hand blowing.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265