

Product Name :
Vernier Calliper

Product Code :
PHYSICSSUK-AS02LAB005



Description :

Vernier Calliper

Technical Specification :

Designed using the best quality raw material, this caliper is used in Physics Lab and Science labs. It is a precision instrument that can be used to measure internal and external distances extremely accurately.

Features:-

- Length: 12.5cm
- Stainless Steel

- Resistant to abrasion
- Accurate measurement
- Robust construction
- Smooth finish.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265