

Product Name :
Aneroid Barometer Demonstration

Product Code :
PHYSICSSUK-AS02LAB006



Description :

Aneroid Barometer Demonstration

Technical Specification :

Aneroid Barometer Demonstration is with 100 mm. dial graduated 28 to 31 inch of mercury and 960 to 1060 millibars. With plated bezel, cover glass and index plastic base approx. 16 cm. diameter. Enclosed in a transparent synthetic case for demonstration purposes. Air tight sealed all round, a pressure bulb is provided to force the movements for educative demonstration.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265