

**Product Name :**  
Balance Double Beam

**Product Code :**  
CHEMISTRYSUK-AS03LAB002C



**Description :**

Balance Double Beam

**Technical Specification :**

Balance Double Beam with single pan and captive weights sliding on scales.

Reading can be taken with index scale and pointer.

Meant for speed in weighing routine samples.

Mounted on wooden base.

Scale Division

On Rear Scale: 0 to 1 gm x 100 mgm.

On Middle Scale: 0 to 1 gm x 100 mgm.

On Front Scale: 0 to 1 gm x 100 mgm.

Capacity: 311 gm.

Sensitivity:  $\pm 2$  mgm.

Distance from suspension: 170 mm.

Hook to pan: 170 mm.

Pan Diameter: 110 mm.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265