

Product Name :
Carboy With Stop Cock

Product Code :
CHEMISTYSUK-AS03LAB004D



Description :

Carboy With Stop Cock

Technical Specification :

These carboys are an excellent replacement for glass carboys and are used for storing and dispensing distilled water & other laboratory solutions.

These polypropylene Carboys are autoclavable, dependable and durable.

These are also equipped with a leak proof spigot screwed to the body.

Material: Polypropylene

Capacity: 10 Ltrs., 20 Ltrs.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265