

**Product Name :**  
Aspirator Bottles

**Product Code :**  
CHEMISTRYSUK-AS03LAB005D



**Description :**

Aspirator Bottles

**Technical Specification :**

Aspirator Bottles are made of Polypropylene and are therefore much lighter than glass. These are autoclavable, durable and dependable. The leakproof spigot screwed to the bottle body requires only a half turn to open or close. These bottles are used for storing or dispensing distilled water & other laboratory solutions. Threads should be totally disengaged before autoclaving.

**Material:** Polypropylene

**Capacity:** 5 Ltrs., 10 Ltrs., 20 Ltrs.



## **Astra Scientific**

[www.astrascientific.com](http://www.astrascientific.com), **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265