

Product Name :
Polypropylene Beaker

Product Code :
CHEMISTRYSUK-AS03LAB001F



Description :

Polypropylene Beaker

Technical Specification :

Beaker, moulded in Polypropylene, have excellent clarity and very good chemical resistance.

Pouring a solution becomes very convenient because of its tapered spout.

These autoclavable Beakers have prominently raised graduations on the outer surface.

Capacity:-

50ml

100ml

250ml

500ml

1000ml

2000ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265