

**Product Name :**  
Wash Bottles New Type

**Product Code :**  
CHEMISTRYSUK-AS03LAB005G



**Description :**

Wash Bottles New Type

**Technical Specification :**

These Bottles are also made of Low Density Polyethylene and are therefore translucent and unbreakable. The distinguishing feature that these bottles contain is the built in delivery tube making it a one piece Wash Bottle. These Bottles are graduated on the outer surface to give an approximate idea about the volume of solution filled or dispensed.

**Material:** Low Density Polyethylene

**Capacity:-**

125ml

250ml

500ml



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265