

**Product Name :**  
Polypropylene Measuring Cylinder

**Product Code :**  
CHEMISTRYSUK-AS03LAB002J



**Description :**

Polypropylene Measuring Cylinder

**Technical Specification :**

Measuring Cylinders, moulded in polypropylene are single piece cylinders.

The board base provides good stability.

These autoclavable cylinders have raised prominent graduations for easy reading and have excellent contact clarity.

**Capacity:-**

10ml

25ml

50ml

100ml

250ml

500ml

1000ml

2000ml



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265