

Product Name :
Measuring Cylinder Hexagonal

Product Code :
CHEMISTRYSUK-AS03LAB003J



Description :

Measuring Cylinder Hexagonal

Technical Specification :

Measuring Cylinders are now available with hexagonal base also.

The hexagonal base provides excellent stability to the cylinder where as the tapered spout works really well while pouring any reagent.

Material: Polypropylene

Capacity:-

10ml

25ml

50ml

100ml

250ml

500ml

1000ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265