

Product Name :
Measuring Cylinder Pentagonal

Product Code :
CHEMISTRYSUK-AS03LAB004J



Description :

Measuring Cylinder Pentagonal

Technical Specification :

Adding to the range of round and hexagonal base cylinders are the pentagonal ones.

Educational Lab Equipments offers you a wide range of cylinders with different bases to choose from.

Prominently raised graduation & excellent contact clarity are among the other features of these cylinders.

Material: Polypropylene

Capacity:-

10ml

25ml

50ml

100ml

250ml

500ml

1000ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265