

**Product Name :**  
Test Tube Stand, Polythene

**Product Code :**  
CHEMISTRYSUK-AS03LAB002M



**Description :**

Test Tube Stand, Polythene

**Technical Specification :**

The base size 210 x 82mm gives stability. Made of polythene high density. Can accept 2 tubes of 25mm dia and four tubes of 16mm dia. 6 vertical pegs at the rear edge of the base ensure quick drying of the tubes. (a) The rack accepts 2 tube of 25mm dia, and 4 tubes of 16mm dia.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265