


Product Name :
Test Tube Stand (Wire Pattern)

Product Code :
CHEMISTRYSUK-AS03LAB004M



Description :

Test Tube Stand (Wire Pattern)

Technical Specification :

These are available in different colors and are made of a special blend of Polypropylene that makes it sink in a water bath and maintain their stability even in agitated water. An economical substitute for wire racks, here are these submersible, auto clavable racks which are stackable when empty. These space saving racks can easily be disassembled for cleaning purposes. The base is in the form of a grid with square openings and has fold up sides and end bars that lock into rack tops.

Material : Polypropylene.

- 13 mm x 90 Tubes.

- 16 mm x 60 Tubes.

- 20 mm x 40 Tubes.
- 25 mm x 24 Tubes.
- 30 mm x 21 Tubes.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265