

Product Name :
Burette Clamps

Product Code :
CHEMISTRYSUK-AS03LAB002P



Description :

Burette Clamps

Technical Specification :

Depressed Burette Clamp made of thick M.S. Sheet, all diepressed symmetrically non-breakable holds and Burette straight two jaws have V-grooves in which Burette rests and other presses it due to spring force.

Well chrome plated.

Bosshead can accommodate upto 18mm dia rod.

For Single Burette.

For Double Burette.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265