

Product Name :
Urine Container

Product Code :
CHEMISTRYSUK-AS03LAB003Q



Description :

Urine Container

Technical Specification :

Moulded in polypropylene these container have excellent contact clarity and the cap has built-in seal that makes it completely leakproof.

These containers can be used for collecting a wide range of samples of sputum, urine etc.

Capacity: 30ml, 60ml.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265