

Product Name :
Ria Vial

Product Code :
CHEMISTRYSUK-AS03LAB005Q



Description :

Ria Vial

Technical Specification :

The polypropylene tubes are clear and can withstand 3000 RCF during Centrifugation while the transparent polystyrene tubes can withstand 1400 RCF.

These polypropylene/ polystyrene tubes are generally used in Ria, Coagulation and Bacteriology.

The PP tubes can accept common acids, alkalies & solvents at room temperature while the PS tubes can accept mild bases, weak acids but not organic solvents.

Capacity:-

12 x 75mm/PP

12 x 75mm/PS

Cap. for 12mm RIA Tubes



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265