

**Product Name :**  
Storage Vial Internal Thread

**Product Code :**  
CHEMISTRYSUK-AS03LAB004R



**Description :**

Storage Vial Internal Thread

**Technical Specification :**

The cap fitted with O-ring screw into the vial makes it a leakproof storage vial.

Another addition to the range of storage vials is this vial which has internal thread in it.

**Capacity:** 2.0ml

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



**Astra Scientific**