

Product Name :
Retort Rings

Product Code :
CHEMISTRYSUK-AS03LAB008M



Description :

Retort Rings

Technical Specification :

A Brass thumbscrew chrome plated to fit boss. Cast iron with boss, hammertone painted. Suitable for fixing on rods up to 13mm dia. The distance between center of ring and supporting rod is approx. 120mm for all sizes.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265