

Product Name :
Test Tube Holder

Product Code :
CHEMISTRYSUK-AS03LAB010M



Description :

Test Tube Holder

Technical Specification :

The jaws will accept tubes up to 38mm diameter and their alignment is maintained by a wire collar. Chrome plated spring steel wire holder with finger grips.

Overall length 130mm.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com