

Product Name :
Float Rack

Product Code :
CHEMISTRYSUK-AS03LAB006P



Description :

Float Rack

Technical Specification :

The surface of the rack has raised alpha numeric index for easy identification and has a central rod convenient lifting.

This Floar Rack moulded in polypropylene is quite compact, autoclavable and can hold 16 micro centrifuge tubes for incubation in water baths.

Easy to load, these are essential in microbiology and immunology labs.

Capacity: 16 MCTs of 1.5ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265