


Product Name :
PCR Tube Rack

Product Code :
CHEMISTRYSUK-AS03LAB008P



Description :

PCR Tube Rack

Technical Specification :

The tube holding spaces are uniformly placed so as to facilitate the use of multipipettes.

Top & base of this rack are made of tough & transparent polycarbonate while the middle portion is made of polypropylene.

Is used for keeping PCR Tubes.

A number of racks can be stacked together to save shelf space.

Capacity: Rack for 0.5ml and 1.5ml MCTs.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265