

Product Name :
Scintillation Vial

Product Code :
CHEMISTRYSUK-AS03LAB006R



Description :

Scintillation Vial

Technical Specification :

These vials can be used with standard counting equipment for liquid scintillation system and are excellent for low activity count.

These unbreakable vials, moulded in polyethylene are fitted with leakproof screw caps.

Capacity: 8ml, 20ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265